



BEHRINGER

Medical gas devices

VAC second stage unit

Product description



The VAC 2nd stage unit is used to shut off VACUUM line pipes; it consists of a Ø 28 mm chemically nickel-plated copper pipe, a 1" ball valve and a G 1/4" F R threaded stub end for the installation of a pressure gauge with electrical contacts and a range of -1/0 bar. Weld-on pipe inlet and outlet connections.

Normatives

UNI EN ISO 9001:2008 | UNI EN ISO 15001:2012

Components

N.One Ø 28 mm chemically nickel-plated copper pipe.

N.One 1" ball valve.

N.One pressure gauge with electrical contacts and a range of -1/0 bar
N.1 G1" 1/4 R - 1" R nipple.

Product codes

GAS	CODE
VAC	HG086

Technical data

CODE	GAS	Ø PIPE
HG086	VAC	Ø28

Behringer S.r.l. via Gualco 4, 16165 Genova ITALY

T +39 010 8309103

F +39 010 8309183

e-mail info@behringer.it

C.F. e P.I. 01167840998